LANDYNE +

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Software	Description of components in the Landyne suite	Reference
PTELS	Periodic table of the elements for the Landyne suite	[2]
SVAT	Structural viewer and analytical tool including atom cluster and layer.	[3]
SPICA	Stereographic projection for interactive crystallographic analysis.	[4]
SAED	Simulation and analysis of electron diffraction (spot) patterns.	[5]
PCED	Simulation of PCED (ring) patterns and phase identification.	[6]
QSAED	Processing, quantification, and analysis of SAED (spot) patterns.	[7]
QPCED	Processing, quantification, and analysis of SAED (ring) patterns.	[8]
HOLZ	Simulation of HOLZ pattern including dynamical correction.	[9]
SMART	Simulation and measurement of rocking curve for crystal thickness.	[10]
SAKI	Simulation and analysis of Kikuchi lines and double diffraction effect.	[11]
TEMUC	Lattice determination of unknown structure in TEM/ED experiments.	[12]
ESPOT ⁺	Electrostatic potential maps derived from electron diffraction patterns.	[10]
CTFscope ⁺	CTF simulation and visualization for conventional and AC-TEM.	[13]
EMIPA ⁺	HRTEM image processing and analysis	[10]
EMCIP ⁺	HRTEM image crystallographic image processing	[10]

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
The dialog box titled "Crystallographic Image Processing" contains a section for "Lattice parameter (nm, °)" with input fields for a = 0.426, b = 0.801, and γ = 90. Below this is a table with four columns: Symmetry, Phase residual (°), R-value (%), and Extinction ratio. The table lists 20 symmetry options, with "p11g" highlighted in blue. At the bottom are buttons for Load, CIP, Save, and Help.

Symmetry	Phase residual (°)	R-value (%)	Extinction ratio
p1	-	-	-
p2	0	-	-
p1m1	11.9	16.1	-
p11m	9.9	16.1	-
p1g1	40.7	25.5	0
p11g	17.1	20.5	0
c1m1			
c11m			
p2mm			
p2mg			
p2gm			
p2gg			
c2mm			
p4			
p4m			
p4g			
p3			
p3m1			
p31m			
p6			
p6m			

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- [14] X.M. Li, H.F. Li, and H.F. Fan. A revised version of the program VEC (visual computing in electron crystallography). Chinese Physics B, 18(6):2459 - 2463, 2009. <http://cryst.iphy.ac.cn/VEC/>.
- [15] Peter Moeck. Crystallographic image processing for scanning probe microscopy. In A. Mendez-Vilas and J. Diaz, editors, Microscopy: Science Technology, Applications and Education, pages 1951 - 1962.
- [16] Jin-Gyu Kim, Seung Jo Yoo, Chang-Yeon Kim, Hyeong-Tae Joul, Circular Fast Fourier Transform Application: A Useful Script for Fast Fourier Transform Data Analysis of High-resolution Transmission Electron Microscopy Image, Applied Microscopy, 2014;44(4):138-143.
- [17] Bart Verberck, Symmetry-Adapted Fourier Series for the Wallpaper Groups, Symmetry 2012, 4, 379-426.

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Copyright 2011-2024 LANDYNE ©